



NATIONAL TYPE EVALUATION PROGRAM

Certificate of Conformance

for Weighing and Measuring Devices

For:

Non-Computing Scale
Electronic Medical
Model: 160-10-7N, 160-10-7NBT, 160-10-7NBLE
 $n_{\max}$: 2750 / 2500
 $e_{\min}$: 0.2 lb / 0.1 kg
Capacity: 550 lb / 250 kg
Platform: 14.4 in x 14.4 in x 1.65 in
Accuracy Class: III

Submitted By: Contact Info. Updated November 2022

Rice Lake Weighing Systems
230 West Coleman St.
Rice Lake, WI 54868
Tel: 715-234-9171
Fax: 715-234-6967
Contact: Brad Bowers
Email: bbowers@ricelake.com
Web site: www.ricelake.com

Standard Features and Options

- Semi-Automatic Zero Setting Mechanism / push button (SAZSM)
- Automatic Zero Tracking (AZT)
- AC to DC power supply converter
- Liquid Crystal Display (LCD)
- Initial Zero Setting Mechanism (IZSM)
- Semi-Automatic Tare (push button)
- DC power supply (6 AA batteries)
- External unit key (kg, lb)
- RS 232 Serial communication port
- USB communication port
- Bluetooth Wireless (optional)
- Low battery indicator/automatic power-off
- Body Mass Index (BMI) is not legal for trade (not NTEP tested)

Load Cell Used:

Model DEBB-300 (Non-NTEP Certified)

Temperature Range: 5 °C to 35 °C (41 °F to 95 °F)

This device was evaluated under the National Type Evaluation Program and was found to comply with the applicable technical requirements of "NIST Handbook 44: Specifications, Tolerances and Other Technical Requirements for Weighing and Measuring Devices." Evaluation results and device characteristics necessary for inspection and use in commerce are on the following pages.

James Cassidy
Chairman, NCWM, Inc.

Kristin Macey
Committee Chair, National Type Evaluation Program Committee
Issued: February 6, 2018

1135 M Street, Suite 110 / Lincoln, Nebraska 68508

The National Conference on Weights and Measures (NCWM) does not approve, recommend or endorse any proprietary product or material, either as a single item or as a class or group. Results shall not be used in advertising or sales promotion to indicate explicit or implicit endorsement of the product or material by the NCWM.



Rice Lake Weighing Systems

Non-Computing Scale / 160-10-7N, 160-10-7NBT, 160-10-7NBLE

Application: Non-Computing digital electronic scale designed for use in the medical field for weight assessment of ambulatory patients.

Identification: The markings are located on a pressure sensitive tamper evident label attached to the back of the indicator housing.

Sealing: The scale is sealed with a category 1 physical seal. Two tamper evident security seals on the left and right side of the indicator housing, affixing a plastic cover over the calibration access opening.

Test Conditions: A Rice Lake Weighing Systems Model 160-10-7BT was submitted for evaluation. The emphasis of the evaluation was on the device design, marking, operation, performance, and compliance with influence factor requirements. For the purpose of the evaluation a 550 lb x 0.2 / 250 kg x 0.1 scale was submitted. Several increasing/decreasing tests, shift tests, and discrimination tests were performed on the scale. Voltage variation tests were conducted at 85 VAC and 264 VAC then at 4.3 VDC and 13.2 VDC. A warm up test, power interruption test, wireless (Bluetooth) test and suitability of level indicator tests were also performed. The scale was tested over a temperature range of 5 °C to 35 °C (41° F to 95 °F). A load of approximately one-half capacity was applied to the scale over 100 000 times. Increase/decrease testing was performed during permanence testing and eccentricity and discrimination testing at the end of permanence testing.

Evaluated By: J. Gibson (OH)



Type Evaluation Criteria Used: *NIST Handbook 44 Specifications, Tolerances, and Other Technical Requirements for Weighing and Measuring Devices*, 2018 Edition. *NCWM Publication 14 Weighing Devices*, 2017 Edition.

Conclusion: The results of the evaluation and information provided by the manufacturer indicate the device complies with applicable requirements.

Information Reviewed By: J. Truex (NCWM)

Examples of Device:

Sealing of Device:



Left Side of indicator

Right Side of indicator